



FENDERWEDGE®



Directa FENDER WEDGE® is designed to protect adjacent teeth during Class II preparation. With a stainless steel protection band, FENDER WEDGE® prevents damage as correct positioning of contact point is established. **Packaged of box of 36 pcs.**

FENDERMATE™



Revolutionary wedge and section matrix combined. Placing of sectional matrix systems and achieving good approximal contacts for Class II restorations can be difficult and challenging. FenderMate offers a unique solution by combining a wedge and sectional matrix in its design.



- Matrix and Wedge inserted as one
- Inserts as easily and rapidly as a wedge
- Right and left available
- Unique wedge design with collapsing wing accommodates most approximal spaces

STERIBAND® RETAINER

Simplifies loading of the matrix band

- Unique clamping device
- Easy loading design
- Eliminates tearing of matrix bands
- Simple dismantling for cleansing
- One retainer for all applications
- Holds all common types of Matrix Bands

Steel Bands are available in 0.04mm width



Luxator® ORIGINAL

**SPECIALLY DESIGNED
SURGICAL CUTTING
INSTRUMENTS**

A real alternative to
traumatic forcep extractions!




Brought to you by:

JS DENTAL
MANUFACTURING, INC.

**Matrix
bands**

Tofflemire Stainless Steel Bands (JR)

144/pkg **\$5.99**



		0.0015"	0.002"
1	Universal	625-1500	625-1540
2	Wide	625-1510	625-1550
3	Narrow	625-1520	625-1560
13	Pedo	625-1530	625-1570

Thinking (Directa)

100/pkg **\$29.95**

Nystrom

Super thin 0.02mm stainless steel

Premolar

X-Narrow	625-1600
Narrow	625-1610
Medium	625-1620

Molar

X-Narrow	625-1630
Narrow	625-1640
Medium	625-1650



Tofflemire (Directa) \$29.95

Super thin 0.02mm stainless steel

1	Universal	625-1670
2	Wide	625-1680
3	Narrow	625-1690
13	Pedo	625-1700



SteriBand, Sigveland, (Directa)

Super thin 0.02mm stainless steel

87mm **\$29.95**


5mm	625-1710
6.5mm	625-1720
8mm	625-1730



T-Bands (Pulpdent)

100/pkg, .002" thick **\$12.99**


	Stainless Steel	Brass
Straight Narrow	625-1900	625-1910
Straight Wide	625-1901	625-1911
Straight Assorted	625-1902	625-1912
Curved Narrow	625-1903	625-1913
Curved Wide	625-1904	625-1914
Curved Assorted	625-1905	625-1915



Ivory Matrix Bands (FCI)

Stainless Steel **\$2.99**
Perforated Bands, 12/pkg

Bicuspid		Molar	
#11	625-0911	#14	625-0914
#12	625-0912	#15	625-0915
#13	625-0913	#16	625-0916



Retainers

Tofflemire

(Arnel) **\$12.99**
Universal Fig 1 630-5000



(First Choice) **\$6.95**
Universal Fig 1 630-1060
Junior Fig 2 630-1061
Senior Fig 3 630-1062



Steriband® Retainer (Directa)

Simplifies loading of the matrix band

630-1070 **\$59.50**



Ivory

Matrix Retainer No 1 (FCI)

Stainless Steel 630-1080 **\$19.50**



Contact Matrix

Contact Matrix (Danville)

Sectional matrix system

Trial Kit 635-2000 **\$49.50**

Contact rings (1 inward & 1 outward)

Thin flex matrices - 50 small & 50 large



Clinical Kit 635-3000 **\$99.50**

Contact rings (3 inward & 3 outward)

Thin flex matrices - 100 small & 100 large

Clinical Kit Plus 635-3010 **\$149.50**

1 Clinical Kit, Ring Plier Plus &

Mega Grip forceps



Ring Plier Plus 635-3020 **\$39.50**

Mega Grip Forceps 635-3030 **\$39.50**

Matrix Placement Tweezers 635-3040 **\$24.50**



Ring Refills - 2 x inward 635-3050 **\$19.50**

Ring Refills - 2 x outward 635-3060 **\$19.50**

MegaRing Refills - 2 x inward 635-3070 **\$29.50**

MegaRing Refills - 2 x outward 635-3080 **\$29.50**

Matrices **\$29.50**

	Thin Flex	Stiff Flex	Ultra Stiff Thin Flex
Small	635-4000	635-4010	635-4020
Large	635-4001	635-4011	635-4021

Wizard Anatomical Wedges (Waterpik)

400/pkg **\$29.50**

Small 640-9160

Medium 640-9161

Large 640-9162

X-Large 640-9163

Assorted 640-9165



Rainbow Wedges (Pulpdent)

Starter Kit, 4 assorted sizes

400/box 640-6200 **\$22.99**

Refill box of 100, all one size **\$6.50**

Small, Blue 640-6210

Medium, Green 640-6211

Large, Yellow 640-6212

X-Large, White 640-6213



Wonder Wedges (Bosworth)

500/pkg **\$29.50**

13mm 640-1000

Asstd 640-1001



Wedges, Barmans (Directa)

Anatomically shaped for better & safer procedure, made from birchwood.

200/pkg **\$15.95**

Mixed 640-3010

Small 640-3011

Large 640-3012



Wedges

Wizard Wedges (Waterpik)

500/pkg **\$29.95**

Assorted 640-9100

Small 640-9120

Large 640-9140

Wide Base 640-9150



Barmans Wedges, Easy Grip (Directa)

Special design provides firm holding & easy safe insertion.

200/pkg **\$15.95**

Large (Easy Grip) 640-3013



Fenderwedges® (JS Dental)

Designed to protect adjacent teeth during class 2 preparation. 36/pkg **\$24.95**

- X-Small, purple 640-3101
- Small, orange 640-3102
- Medium, green 640-3103
- Large, yellow 640-3104
- Assorted, 144/pkg 640-3100 **\$79.50**



Fendermate™ (JS Dental)

Revolutionary wedge & section matrix combined.



- Assorted, 72/pkg 640-3200 **\$149.50**
- Individual 18/pkg **\$39.95**
- Left Reg 640-3201
- Left Narrow 640-3202
- Right Reg 640-3203
- Right Narrow 640-3204
- Clamping Ring
- Short 640-3205
- Long 640-3206

Wykle Wedges (Wykle Research)

- 500/box **\$17.50**
- 11mm Natural 640-8000
- 13mm Yellow 640-8001
- 15mm Blue 640-8002
- 17mm Green 640-8003
- Assorted 640-8004



Plastic Wedges (Almore)


- 250/pkg **\$15.95**

Uni-Wedge

- Regular Orange 640-2000
- Narrow Green 640-2001

Bi-Wedge

- Regular Orange 640-2010
- Narrow Green 640-2011



Fintrec Transparent Wedges (Pulpdent)

Silver coated, advanced wedge design, tapers on 4 sides so that one size fits all, .05mm - 2mm

- 250/pkg 640-5000 **\$33.95**



Contact Wedges (Danville)

Elastic, create ideal cervical seal, fast & easy to use, eliminate hash & ledging

- Kit - Small, Medium, Large 255 pcs 640-1500 **\$74.50**
- Refill 85/pack **\$29.50**
- Small 640-1501
- Medium 640-1502
- Large 640-1503



Matrix Strips

Mylar Matrix

Plastic for composite

- 50/pkg 640-1600 **\$1.99**
- 100/pkg 640-1610 **\$2.99**



Directa Strips (Directa)

Dispenses one strip at a time

- 100/pkg **\$6.99**
- Straight 640-1620
- Curved 640-1630



Root Form (Directa)

Cervical Matrix
Simplifies filling of excavated subgingivally exposed roots

- Kit 640-1640 **\$34.95**
15 x 6mm & 15 x 8.5mm
- Refill 640-1641 **\$19.95**
15 x 6mm
- Refill 640-1642 **\$19.95**
15 x 8.5mm



Matrix Strips (Crosstex)

Plastic for composite

- 500/pkg **\$4.99**
- Long 640-1700



* Prices are subject to change without notice

